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## INSTRUCTION FOR USE

### 1. Name

NITHYA – Powder of equine collagen Type I, sterile.

### 2. Composition of medical device and accessories

The medical device is composed at 100% of collagen powder, sterile; pyrogen-free. The limit content of bacterial endotoxins not more than 0.5 U / ml.

The CAS number of collagen powder is: 9007-34-5.

Type I collagen is characterized by its complete and complex molecular form: a protein in the typical triple-helical structure.

### 3. Intended use

Aesthetic surgery. Chrono and photo aging of the face, neck and hands.

### 4. Packing set

Each pack contains n.3 of glass vials each containing 70 mg of equine Type I collagen powder, sterile.

The packaging provides 1 leaflet, too.

The needles are not provided in package, but 30-32 G ½" – 4 mm needles are recommended. Syringes are not provided in package; 5 syringes of 1 mL are recommended.

The choice of syringes and needles is at discretion of physician.

### Packaging materials:

| No | Item                      | Material | Composition                                                                                                                                                                          | Grade                   | Size                                                                                                                                                     | Supplier                                                                                           |
|----|---------------------------|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| 1. | Vial<br>capacity<br>10 ml | Glass    | SiO <sub>2</sub> 71-73%<br>Al <sub>2</sub> O <sub>3</sub> 2-4%<br>CaO 10-12%<br>K <sub>2</sub> O 0-2%<br>Na <sub>2</sub> O 12-14%<br>Fe <sub>2</sub> O <sub>3</sub> 0-1%<br>MgO 0-1% | Yellow glass<br>type 3. | height 53 mm ± 1 mm;<br>the diameter of body 25 mm ± 1 mm;<br>outer diameter of the neck 20 mm ± 1.0 mm;<br>inner diameter of the neck 13.0 mm ± 0.5 mm; | SIMONETTA<br>SAS.<br><br>Via G. Tartini,<br>8 - 20158,<br>Milan - Italy                            |
| 2. | Bung                      | Rubber   | Chlorobutyl                                                                                                                                                                          | PA1938                  | cap diameter 20 mm ± 1 mm;<br>base diameter 13.0 mm ± 0.5 mm;<br>height 9.0 mm ± 0.5 mm;                                                                 | PHARMA<br>TRADE S.R.L.<br><br>Zona Ind.le –<br>Loc. Punta<br>Penna- 66054<br>VASTO (CH) -<br>Italy |

|    |              |           |                 |        |                                                                                                                                    |                                                                                                    |
|----|--------------|-----------|-----------------|--------|------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| 3. | Flip-off cap |           |                 |        |                                                                                                                                    | PHARMA<br>TRADE S.R.L.<br><br>Zona Ind.le –<br>Loc. Punta<br>Penna- 66054<br>VASTO (CH) -<br>Italy |
|    | - seal       | Aluminium | Aluminium alloy | 8011   | outer diameter<br>(assembly) 23.0 mm $\pm$ 0,5 mm;<br>height (assembly)<br>8.0 mm $\pm$ 1,0 mm;<br>diameter<br>20 mm $\pm$ 0,5 mm; |                                                                                                    |
|    | -cap         | Plastic   | Polypropylene   | HP456H | outer diameter 23.0 mm $\pm$ 0,5 mm;<br>inner diameter 20,0 mm $\pm$ 0,5 mm;<br>height 3.0 mm $\pm$ 0,3 mm;                        |                                                                                                    |

## 5. Technical characteristics

- Weight of the powder 70 mg  $\pm$  5 mg;
- Total weight of finished vial that: 25g  $\pm$  10%;

## 6. Risk class

Class 3.

## 7. Indications

Chrono and photo aging of the face, neck and hands. The use to the body can be extended according to the medical discretion. Indicated for the regeneration of the dermal and subdermal connective, through a mechanism of fibroblastic stimulation and neocollagenesis. Remodeling and genetic malformation such as, for example: hypotrophy and skin hypotonia, reduced skin elasticity.

## 8. Contraindications

- Do not inject NITHYA in the periorbital region (eyelid, crow's feet) and glabellar (front).
- The use of NITHYA in the area of the circles is reserved for physicians who have received training specific to this technique and possess a good knowledge of the physiology of this region.
- Do not inject into the blood vessels (intravascular).
- Do not exceed corrective action.

It should not be used:

- In patients with epilepsy not controlled by treatment;
- In patients with a tendency to develop hypertrophic scars;
- On patients with a known hypersensitivity to collagen;
- In patients with porphyria;
- During pregnancy or breast-feeding;
- On children;

- Areas of skin inflammatory and / or infectious (acne, herpes, etc.);
- In immediate association with a laser treatment, a deep chemical peel or dermabrasion. In the case of cosmetic peeling, it is recommended not to inject if the inflammatory reaction is important.

## 9. Possible side effects

The product is hypoallergenic. There are no reported side effects.

The patient must be informed of the existence of secondary potential effects related to this device that may appear immediately or after a certain period of time.

Among them, there are reported (non-exhaustive list):

- Inflammatory reactions (redness, swelling, erythema, etc.), sometimes associated with itching and pain to the touch, or both, which may occur after the injection. Such reactions may persist for a week.
- Haematoma.
- Hardening or nodules at the point where the injection has been practiced.
- Staining or discoloration of the injection area.
- Poor treatment efficacy or poor filling effect.
- Cases of necrosis in the glabellar region abscesse, granulomas and immediate or delayed hypersensitivity have been described in the literature after collagen injection. It is therefore necessary to take into account these potential risks.
- The persistence of inflammatory reactions over the week or the appearance of other side effects must be promptly reported by the patient at the doctor, who will be required to eliminate these reactions through a suitable care.
- Any other unwanted secondary effect related to the injection of "product" should be reported to the distributor and/or the manufacturer.

## 10. Mode of application

Dilute the collagen powder contained in the vial from 70 mg to obtain a suspension in 5 ml of physiological solution or distilled water.

NITHYA has to be brought in suspension immediately before use and should be considered single use. Administered with small gauge needles (30G) at the level of the deep dermis or dermal junction hypodermic, through micro wheals. Distribute the product in the region to be treated and then run a light massage.

Surgical Technique: the medical device is exclusively intended for insertion by skilled surgeons.

## 11. Operating procedure

Administered with small gauge needles (30G) at the level of the deep dermis or dermal junction hypodermic, through micro wheals.

Surgical Technique: the medical device is exclusively intended for insertion by skilled surgeons.

## 12. Indications for use

Inject slowly. Failure to follow these instructions can result in the risk of needle leakage and/or loss of product in the Luer-lock level. The amount of product to be injected will depend on the type of area to be corrected. After the injection, it is important to massage the area treated so as to ensure a uniform distribution of the product.

### 13. Precautions of use

Use normal precautions for intradermal infiltration.

As with any percutaneous treatment using the product may present a risk of infection.

Treatment can result in a slight stinging which disappears in a few minutes. The treated area may have reactions commonly associated with injecting trauma, such as swelling and erythema with spontaneous resolution within a few hours after treatment.

Do not use in case of known hypersensitivity to the product.

The treatment is not indicated in patients with collagen vascular and autoimmune disease of the connective tissue (connective tissue).

Do not inject into areas where the skin infection or inflammation are present in place.

Use during pregnancy and lactation is not recommended, although there are no known specific contraindications.

Surgical Technique: the medical device is exclusively intended for insertion by skilled surgeons.

The doctor's technical skills is essential for the success of the treatment, and the device should be used by physicians who have received training specific to the filling injections technique.

Any modification or use of this product, other than the conditions of use defined herein, may impair its sterility and its performance, which can no longer be guaranteed. Before undertaking the treatment, it is recommended to inform the patient of the indications and contraindications of the device itself, its incompatibilities and the side effects that could occur.

### 14. Methods of cleaning

Before injection, proceed to thoroughly disinfection the area to be treated.

### 15. Sterility

Gamma irradiation.

### 16. Storage conditions and Terms of use

To store NITHYA at a range temperature between 0°C and +30°C, in dry environment and in its original packaging. Protect from sunlight, heat and handle with care.

Humidity not more 60%. Transportation at range of temperature between 0°C and +30°C, 60% humidity.

### LIST OF APPLICABLE STANDARDS

The medical device "Collagen powder vials, sterile" having the brand name NITHYA obtained the EC Certificate in compliance with the following national and international applicable standards:

- 93/42/EC Directive concerning the medical device, recognised in Italy by D.Lgs n. 46 of 24.02.1997 e ss.mm.ii;
- 2007/47/EC Directive of 05.09.2007 amending the 93/42/EC Directive, recognized in Italy by D.Lgs n. 37 of 25.01.2010;
- U NI EN ISO 14971/2012 - Application of risk management to medical devices; ICH Q9/2011-Quality risk management;
- UNI EN ISO 22442/2008 - Medical devices utilizing animal tissues and their derivatives;
- U NI CEI EN 1041/2009 - Information supplied by the manufacturer of medical devices;
- UNI CEI EN 15223-1:2012 - Symbols to be used with medical device labels, labelling and information to be supplied;
- U NI EN ISO 13485/2012 - Quality management systems - Requirements for regulatory purposes;
- UNI EN ISO 10993-5:2009 and U NI EN ISO 10993-10:2010- Biological evaluation of medical devices - Biocompatibility tests;
- U NI EN ISO 11137:1, 2, 3 - Sterilization of Medical Devices.

**17. Shelf Life**

2 years.

**18. Labelling****Consumer packaging:**

| NAME                                                                                | NITHYA                                                                              |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Composition                                                                         | Type I Collagen                                                                     |
| Manufacturer                                                                        | Euroresearch S.r.l. – Via Larga 15 – 20122 Milano<br>– Italy;<br>Tel+ 39 02 8055660 |
|    | Do not reuse!                                                                       |
|    | Do not use if packaging is damaged!                                                 |
|   | Reading information contained within the box                                        |
|  | Storage temperature from 0 °C to 30 °C                                              |
|  | Sterile Product: sterilized by irradiation                                          |
|  | Product Expiry Date                                                                 |
| CE xxxx                                                                             | CE Mark                                                                             |
| xxxx                                                                                | LOT number                                                                          |
| month/year                                                                          | Shelf life                                                                          |

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**19. UTILIZATION OF MEDICAL DEVICE****Utilization**

Products that have expired should be destroyed as medical waste of Class A (epidemiologically non hazardous waste, close in composition to the municipal solid waste).

Disposal of packaging materials should also be carried out in accordance with standard practice in the medical institution, and in the manner prescribed by local, federal and state rules and regulations.

**20. WARRANTY**

The manufacturer, Euroresearch S.r.l. (from Italy), ensures the compliance of medical device with the requirements of the Technical File, if to follow conditions of use, storage and transportation. The guaranteed shelf life is 2 years from the date of manufacturing. The product in the manufacturer's packaging should be stored in the warehouses of the supplier at a storage temperature of 0 ° C to 30 ° C in a dry place protected from light, humidity no more than 60%. Manufacturer assumes absolute responsibility in case of defective goods and replace them.

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With regard to the recorded request,  
submitted by the legal representative/proxy  
of the Enterprise, we declare that above mentioned  
person has the power to sign documents  
for international trade exchanges

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This Apostille only certifies the signature, the  
capacity of the signer and the seal or stamp i  
bears. It does not certify the content of the  
document for which it was issued

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